

FAM108A1-family with sequence similarity 108, member A1 Gene

Reference: PTXBC000158

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	family with sequence similarity 108, member A1
GENE ID	81926
ORIGIN SPECIES	Human
ACCESSIONS	BC000158
VECTOR	pDONR223
NCBI SYMBOL	FAM108A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatggcgtgctgagtgagctctgctcctctctgctcccgccctgcccggcgcatcgctccaagctcgctctctgcccggaggccactactcctgggtgctgagcccgagccgggctgtggggccggggcgccccctggggaccctgagagcctct
cggcgccacccggcgctggaagctgcacctgacggagcgtgcccgaattccagtagccagcgagctggacaccatcgaggtctcccccacaaagagcccgccggcaaccagctctctgcatgtatgttcgctgctgctgctggccaggtacacggtctctctcgacgg
caatgccgtggacctgggcccagatgagcagctctacattggcctgggctcccgctccactgcaacatctctctacgactactcgggtacgggtgccagctcggcagggccttcgagaggaacctctatgccgacatcgaccccgctggcagggcctgcgcaccaggtacggca
tcagcccgacagcatcctgtacgggagagcatggcacgggtgcccacgtggactggcctcgctacgagtggtgctgctgcactcgccgctcactcgggcatgcgctgcttccccgacaccaagaagacctactgttcgacgcttccctaagtga