

## KRT18-keratin 18 Gene

Reference: PTXBC000180

|                |                                                                                                                                 |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CK-18; CYK18; K18; keratin, type I cytoskeletal 18; cell proliferation-inducing gene 46 protein; keratin 18, type I; keratin 18 |                     |                |

**For research use only.**

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | keratin 18                                                                                       |
| <b>GENE ID</b>          | 3875                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC000180                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | KRT18                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgagcttcaccactcgtccacctctccaccaactaccggtccctgggtctgtccaggccccagctacggcgccggccggtcagcagcgccgccagcgtctatgcaggcgctgggggctctgtgtcccgatctccgtgtcccgctccaccagcttcaggggcgccatggggtc  
cgggggcttgccaccgggatagccgggggtctggcagggaatgggggcatccagaacgagaaggagaccatgcaaaagcctgaacgaccgctggcctcttacctggacagagtggagagcctggagaccgagaaacggagctggagagcaaaatccgggagcacttggag  
aagaagggacccaggtcagagactggagccattctcaagatcatcgaggacgtgagggtctcagatcttcgcaatactgtggacaatgccgcacgtctctgcagattgacaatgccgtcttctgctgatgactttagagtaagatgagacagagctggccatgcgccagtctgt  
ggagaacgacatccatgggtctcgcgaaggtcattgatgacaccaatatacacgactgcagctggagacagagatcgagggtctcgaaggaggagctgctctcatgaagaagaaccacgaagaggaaagtaaaaggcctacaagccagattgccagctctgggttgaccgtggaggtgta  
gatgcccccaaatctcaggacctcgccaagatcatggcagacatccgggccaatatacagagctgggtcgggaagaaccgagagagctagacaagtactggtctcagcagattgaggagagcaccacagtggtcaccacacagctgctgaggttggagctgctgagacgacgctc  
acagagctgagagctacagctccagtccttggagatcgacctggactccatgagaatctgaaggccagcttggagaacagcctgaggagggtggaggccgctacgcccctacagatggagcagctcaacgggatcctgctgcaccttgatcagagctggcacagaccgggcaga  
gggacagcgccaggccaggaatgatgaggccctgctgaacatcaaggtcaagctggaggctgagatcgccacctaccgcccgtctggaagatggcgaggacttaattcttggtgatgccttggacagcagcaactcatgcaaacatccaaaagaccaccaccgcccgatagt  
ggatggcaagtggtgctgagaccaatgacacaaagtctgagcattaa