



Reference: PTXBC000185

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CPT1; CPT1-L; L-CPT1; carnitine O-palmitoyltransferase 1, liver isoform; CPT I; CPTI-L; carnitine O-palmitoyltransferase I, liver isoform; carnitine palmitoyltransferase 1A (liver); carnitine palmitoyltransferase I, liver; carnitine palmitoyltransferase 1		

For research use only.

GENE DESCRIPTION	carnitine palmitoyltransferase 1A (liver)
GENE ID	1374
ORIGIN SPECIES	Human
ACCESSIONS	BC000185
VECTOR	pDONR223
NCBI SYMBOL	CPT1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]