

MRPL51-mitochondrial ribosomal protein L51 Gene

Reference: PTXBC000191

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDA09; HSPC241; MRP64; bMRP64; 39S ribosomal protein L51, mitochondrial; L51mt; MRP-L51; bMRP-64; mitochondrial ribosomal protein 64; mitochondrial ribosomal protein bMRP64; mitochondrial ribosomal protein L51		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L51
GENE ID	51258
ORIGIN SPECIES	Human
ACCESSIONS	BC000191
VECTOR	pDONR223
NCBI SYMBOL	MRPL51
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagggaacctctatccgggaggttaggcgcctgtgggactgggtgcctctggcgtgcagaagcttctctgtgtgctagattgatcggtataaggctcactctccgcccccaaaagtgggtgatcgttgaacgagaaaaggccatgttcggagtgtatgacaacatcgg
gatcctgggaaactttgaaaagcaccaccaagaactgatcaggggcccatatggcttcgaggttgaaaagggaatgaattgcaacgttgtatccgaaagaggaaaatggttggaagtagaatgttcgtgatgacctgcacaaccttaataaacgcacccgctatctctacaacacittaa
ccgacatgggaagtctcgtatag