

CCNG1-cyclin G1 Gene

Reference: PTXBC000196

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCNG; cyclin-G1; cyclin-G; cyclin G1		

For research use only.

GENE DESCRIPTION	cyclin G1
GENE ID	900
ORIGIN SPECIES	Human
ACCESSIONS	BC000196
VECTOR	pDONR223
NCBI SYMBOL	CCNG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgataagggtactgacaacaactgactctcagaaactgctacaccagctgaatgccctgttggacaggagctctagatgtcagccaaaggtctgtgttggactaattgagctctgcacacgataatggcctcagaatgactgcaagactaagggtcttgaagtaaaagatctttagtct
aactcagttcttggcttgacacagagacattttcttagctgtgaatttactggacagatctctgtctaaaatgaagggtacagcccaagcacttgggtgtgttggactgagctgttttatttggctgtataaatcaatagaaggaaaggaatgtccattggcactgactgatccgaataag
tcaatataggtttacggtttcagacttgatgagaatggaagaattgtattggagaaggtgtgttggaaagtcaaaagctactgctttcaattctgcaactgtatttactcctcagagaaacttgcacttgaagggaagaatagcattatttgaagactagaagctcaactgaaggc
atgtcattgcaggatcatattttctaaagcaaaccttctgtgttggcattgtctatcattgcattagatccaagcacagaagtgtagagtaacagaaggaatagaatgtcttcagaaacattccaagataaatggcagagatctgaccttctggcaagcctgtatccaaatgtttaaactga
atattcatcaataaagtggtccaaaccaaatgttcagaagtgaaatggattgtttctggcgtactgcacggcaattgaagcatagctactacagaataactcacctccaacaattctgaaatggctccttaa