

RPS3A-ribosomal protein S3A Gene

Reference: PTXBC000204

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FTE1; MFTL; S3A; 40S ribosomal protein S3a; fte-1; v-fos transformation effector protein 1; ribosomal protein S3A		

For research use only.

GENE DESCRIPTION	ribosomal protein S3A
GENE ID	6189
ORIGIN SPECIES	Human
ACCESSIONS	BC000204
VECTOR	pDONR223
NCBI SYMBOL	RPS3A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggtggcaagaacagcgccctacgaaagcgggcaaaaaggagcgaagaagaagtgggtgacatttctaaagaaagattggatgatgtaaacacctgctatgttcaataaagaattatggaaagacgctcgtcaccaggaccaagggaacaaaatgcacatgatg
gtctcaagggtcgtgtgtttgaagtgaagctctgctgatttgcaaatgatgaagtgcatttagaaaaatcaagctgattactgaagatgttcagggtaaaaactgctgactaacctccatggcatggatctaccggtgacaaaatgtgtccatgtcaaaaatggcagacaatgattgaagct
cacgttgatgtcaagactaccgatggttactgtctctgtctgtgtgtgttttactaaaaaacgcaacaatcagatcacggaagacctctatgctcagcaccaacaggtccgcaaatccggaagaagatgatggaatcatgacccgagaggtgcagacaatgacttgaagaagtgg
tcaataaattgattccagacagcattggaaaagacatagaaaaggcttccaatctattatcctcctcatgatgtctcgttagaaaagtaaaatgctgaagaagcccaagttgaattgggaaagctcatggagcttcattggtgaaggcagtagttctgaaaagccactggggacgagaca
gggtcctaaagttgaacgagctgatggatatgaaccaccagtccaagaatctgttaa