

THAP4-THAP domain containing 4 Gene

Reference: PTXBC000247

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-36; PP238; THAP domain-containing protein 4; THAP domain containing 4		
For research use only.			

GENE DESCRIPTION	THAP domain containing 4
GENE ID	51078
ORIGIN SPECIES	Human
ACCESSIONS	BC000247
VECTOR	pDONR223
NCBI SYMBOL	THAP4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccatcaacgaggtcatcctgtcggcgtcagggcctgcaagctcatcgaactcactgactcctctctcctccggcagaacaagagccaggtgtgtcctcgtcgggagcaggtggagaagaagaacgagcgtgaagacctgcggcagagggtcagccgctccgac
agccaggtgcgggaagctacaggagaagctggatgagctgaggagagtgagcgtccctatccaagtagcctgctgtcggccagccgagcccccaagatgaaccagtggtggagccactgtcctggatgctgggcacctgtgtcggaccacctggagccgggacctaccc
cacactgcagcccttcagctacgtgaggaggttcacatctccacgtggccagcccatgctgaactctcgttcaactcctccacccggacacgcgaagccgatgcacagagagtggtgcttcattcgcctcaagccgacaccaacaaggtggccttgcagcggccagaacac
aggcgtggtggaagtgaggaggggcaggtgaacggcaggagctgtcatcgcacccatccatgccaggatctcttgcgaaggagccccacgtagagcagatcacccggaaggtcaggctgaattctgaaggcaactgagcagacggtctccatggcaaccacgacac
agccaatgactcagcatcttcagctcactacaagaagtgaccccgtaa