



Reference: PTXBC000249

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CSNB3; CSNBAD2; PDEB; RP40; rd1; rod cGMP-specific 3',5'-cyclic phosphodiesterase subunit beta; GMP-PDE beta; phosphodiesterase 6B, cGMP-specific, rod, beta; rod cGMP-phosphodiesterase beta-subunit; phosphodiesterase 6B		

For research use only.

GENE DESCRIPTION	phosphodiesterase 6B, cGMP-specific, rod, beta
GENE ID	5158
ORIGIN SPECIES	Human
ACCESSIONS	BC000249
VECTOR	pDONR223
NCBI SYMBOL	PDE6B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr