



Reference: PTXBC000254

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ACTRIB; ACVRLK4; ALK4; SKR2; activin receptor type-1B; activin A receptor, type IB; activin A receptor, type II-like kinase 4; activin receptor-like kinase 4; serine/threonine-protein kinase receptor R2; activin A receptor type 1B		

For research use only.

GENE DESCRIPTION	activin A receptor, type IB
GENE ID	91
ORIGIN SPECIES	Human
ACCESSIONS	BC000254
VECTOR	pDONR223
NCBI SYMBOL	ACVR1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr