

CUEDC2-CUE domain containing 2 Gene

Reference: PTXBC000262

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C10orf66; bA18I14.5; CUE domain-containing protein 2; CUE domain containing 2		
For research use only.			

GENE DESCRIPTION	CUE domain containing 2
GENE ID	79004
ORIGIN SPECIES	Human
ACCESSIONS	BC000262
VECTOR	pDONR223
NCBI SYMBOL	CUEDC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagctggagagatcgtcagtcagccctccttgccttgcagacacacctcccgaggccgacctcagtcgcttgatgaggtcatcttctctatgtgcttgggtcctggaggacctgggcccctcgggcccacagaggagaactcgatatggagcttctactgagatgatggagccctatgtgcttgcctcggccacatccccaggggcacaaataggggacatgatgcagaagctctcaggcagctgagcgaatgccaggaacaaagagaacctgcaaccgcagagctctgttccaaagtcaggtgcccctctcccagagccccctgcagcggccgaaatgctcaaagaagagactaggtcttcgctgctgctgcagacaccaagatgaggcaactggcgtgaggaggagctctgccagggtggatgtactctggaggtgttccctacctgttcggtggagcagcccagtggtgctggccaaagctcgggggactggagaagctgtgcagatgctgtagagggaaggaaggggcctgcagcctggaggggcccaaccaggaacctgccagagcctcagaggccccaaaaggatgagctgaagtccttcctcgcagaagtacatgatggtgatagcgcatag