

RNF220-ring finger protein 220 Gene

Reference: PTXBC000279

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	E3 ubiquitin-protein ligase RNF220; C1orf164; ring finger protein 220		
For research use only.			

GENE DESCRIPTION	ring finger protein 220
GENE ID	55182
ORIGIN SPECIES	Human
ACCESSIONS	BC000279
VECTOR	pDONR223
NCBI SYMBOL	RNF220
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttacaatgtgtgtataaatgggacaactctcgcctctacctgtccctcccttgggtgtatgattttcttttttaagaacctggagcagtgccctctcagggttggtggagctcgcccatccacctcttgggtgtatctgcctctctctctctgtgtgtccctccctctcc
catgtctcgggtgtcagtggtgtatattctctccagacatggggcacacgcccaaggacatgatcctctccttagctctagctcatgggctctttataaggagtgggggtagagagcaggaaatgggaaccgagctga