

ARD1A-ARD1 homolog A, N-acetyltransferase (*S. cerevisiae*) Gene

Reference: PTXBC000308

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARD1A; ARD1; ARD1P; DXS707; MCOPS1; NATD; OGDNS; TE2; N-alpha-acetyltransferase 10; ARD1 homolog A, N-acetyltransferase; N-acetyltransferase ARD1, human homolog of; N-terminal acetyltransferase complex ARD1 subunit homolog A; natA catalytic subunit Naa10;		

For research use only.

GENE DESCRIPTION	ARD1 homolog A, N-acetyltransferase (<i>S. cerevisiae</i>)
GENE ID	8260
ORIGIN SPECIES	Human
ACCESSIONS	BC000308
VECTOR	pDONR223
NCBI SYMBOL	ARD1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacatccgcaatgcagagccagaggacctaataacatgcagcactgcaacctcctctgctgcccagaaactaccagatgaataactacttaccatggccttctgcccagctctctacattgctgaggacgagaatgggaagattgtgggtatgtcctggccaaatgga
agaggacccagatgatgtgcccacggacatatcacctattggtgtgaagcgttcccacggcgctcgtctggtcagaaactgatggaccaggcctctcgagccatgatagagaactcaatgccaaatattgtccctgcatgcaagagtaaccggcgccctgcacctc
tatccaacacctcaacttcagatcagtgaaatggagcccaatactatgcagatgggagagacgcctatgccatgaagcgggacctcactcagatggccgacgagctgaggcggcacctggagctgaaagagaaggcagggcagctggtgctgggtgccatcagaacaagggtg
gagagcaaggcaattcacctccgagctcagagagagcctgtcgcgaggagaaggcctggtgcccagagatgtgtgggacagcaaggacctcagcgaggtcagcgagaccacagagacagatgtcaaggacagctcagaggcctccgactcagcctctag