



Reference: PTXBC000326

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	beta'-COP; coatomer subunit beta'; beta'-coat protein; betaprime-COP; coatomer binding complex, beta prime subunit; coatomer protein complex subunit beta prime; coatomer protein complex, subunit beta 2 (beta prime); p102; coatomer protein complex subunit		

For research use only.

GENE DESCRIPTION	coatomer protein complex, subunit beta 2 (beta prime)
GENE ID	9276
ORIGIN SPECIES	Human
ACCESSIONS	BC000326
VECTOR	pDONR223
NCBI SYMBOL	COPB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]