



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

FADD-Fas (TNFRSF6)-associated via death domain Gene

Reference: PTXBC000334

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

GIG3; MORT1; FAS-associated death domain protein; Fas (TNFRSF6)-associated via death domain; Fas-associating death domain-containing protein; Fas-associating protein with death domain; growth-inhibiting gene 3 protein; mediator of receptor-induced toxicit

For research use only.

GENE DESCRIPTION	Fas (TNFRSF6)-associated via death domain
GENE ID	8772
ORIGIN SPECIES	Human
ACCESSIONS	BC000334
VECTOR	pDONR223
NCBI SYMBOL	FADD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr