

RPS28-ribosomal protein S28 Gene

Reference: PTXBC000354

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DBA15; S28; 40S ribosomal protein S28; ribosomal protein S28		

For research use only.

GENE DESCRIPTION	ribosomal protein S28
GENE ID	6234
ORIGIN SPECIES	Human
ACCESSIONS	BC000354
VECTOR	pDONR223
NCBI SYMBOL	RPS28
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacaccagccgtgtgcagcctatcaagctggccaggtcaccaaggtctctgggcaggaccggttctcaggggacagtgacgcaggtgcgcgtggaattcatggacgacagcagccgatccatccgaatgtaaaagggcccggtgcgcgagggcgacgtgtcaccccttttg
gagtcagagcagagaagcccgaggttgcgtga