



Reference: PTXBC000392

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	POR1; arfaptn-2; partner of RAC1 (arfaptn 2); ADP ribosylation factor interacting protein 2		

For research use only.

GENE DESCRIPTION	ADP-ribosylation factor interacting protein 2
GENE ID	23647
ORIGIN SPECIES	Human
ACCESSIONS	BC000392
VECTOR	pDONR223
NCBI SYMBOL	ARFIP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1