

CCDC106-coiled-coil domain containing 106 Gene

Reference: PTXBC000412

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HSU79303; ZNF581; coiled-coil domain-containing protein 106; protein predicted by clone 23882; coiled-coil domain containing 106		

For research use only.

GENE DESCRIPTION	coiled-coil domain containing 106
GENE ID	29903
ORIGIN SPECIES	Human
ACCESSIONS	BC000412
VECTOR	pDONR223
NCBI SYMBOL	CCDC106
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaatgacgggagcagtcggaggcgacaatgaaggacgatgagaccttcgagatcctcccttcgatgaggcaccacacagatctttacagtcctgagccctctcgagaaactcgaggagcctcggaggctgcgtcctccgctgctgatgaacagcgtcaagacccagctgcacatggctctggagagggaactcctggctgcagaagcgcatcgaggacctggaggagagaggacacttcctgctgctggatggaggcagaggaccactgccggatgaagcctgggccaggcggaatggaggggagacagcctgggtgggctggggcgaggcctcggaacctgagtcagcagcctcctccctcagcggagcgtccgaagaaggcagtgccagtgagaggaggcgccagaagcagaaggagggtgctagtcgaggcgctttgggaagcccaaggccgggagaggcagcgagtgaaggacggcgacgggtcctctgcccgtacaagaagatcctgggcaccttcagaagctcaagagcatgtcggggccttcgagcaccacgcgtggacaggaacacgtggcgtgaccacgcccacgcggagctgctcattgtggccccgagaagctggccgaggtggcgagttcgaccctccaaggagcgctgctgagtaactccgctgctttctggccctggagcagagacgctcaagaaggctcagggcgctcaagaagagcaagctgctgctgccatcacctaccgcttcaagcggtga