



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C1QBP-complement component 1, q subcomponent binding protein Gene

Reference: PTXBC000435

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

GC1QBP; HABP1; SF2p32; gC1Q-R; gC1qR; p32; complement component 1 Q subcomponent-binding protein, mitochondrial; ASF/SF2-associated protein p32; C1q globular domain-binding protein; complement component 1, q subcomponent binding protein; glycoprotein gC1q

For research use only.

GENE DESCRIPTION	complement component 1, q subcomponent binding protein
GENE ID	708
ORIGIN SPECIES	Human
ACCESSIONS	BC000435
VECTOR	pDONR223
NCBI SYMBOL	C1QBP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr