

MIF-macrophage migration inhibitory factor (glycosylation-inhibiting factor) Gene

Reference: PTXBC000447

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GIF; GLIF; MMIF; macrophage migration inhibitory factor; L-dopachrome isomerase; L-dopachrome tautomerase; phenylpyruvate tautomerase; macrophage migration inhibitory factor (glycosylation-inhibiting factor)		

For research use only.

GENE DESCRIPTION	macrophage migration inhibitory factor (glycosylation-inhibiting factor)
GENE ID	4282
ORIGIN SPECIES	Human
ACCESSIONS	BC000447
VECTOR	pDONR223
NCBI SYMBOL	MIF
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgatgttcacgttaaacaccaacgtgccccgcctccgtgccggagcgggttctctccgagctaccacagcagctggcgaggccacgcggcaagccccccagtcacatcggtgcacgtggtccggaccagctatggccttcggcggtccagcagccgtgcgcgtc
tgcagcctgcacagcatcggcaagatcggcgcgcgagaaccgctctacagcaagctgtgtgcggcctgctggccgagcgctgcgcacatcagccggagaggggtctacatcaactattacgacatgaacgcggccaatgtgggtggaacaactccacctgcctaa