

TSC22D1-TSC22 domain family, member 1 Gene

Reference: PTXBC000456

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Ptg-2; TGFB1I4; TSC22; TSC22 domain family protein 1; TGFB-stimulated clone 22 homolog; TGFbeta-stimulated clone 22; cerebral protein 2; regulatory protein TSC-22; transcriptional regulator TSC-22; transforming growth factor beta-1-induced transcript 4 pr		

For research use only.

GENE DESCRIPTION	TSC22 domain family, member 1
GENE ID	8848
ORIGIN SPECIES	Human
ACCESSIONS	BC000456
VECTOR	pDONR223
NCBI SYMBOL	TSC22D1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaatcccaatgggtgtagaccagtgccgatgtaggtttaccaactgagacattttcaatttcttctgtcatccttctgtgggactgaaaacgttctgtgagacttgataatagctccttgggtgcaagtgtggtgtagctattgacacaaatcgagcaagctatggatctagtgaagccatttgatgtatgcggtcagagaagaagtggaggtcctcaagagcaaatcaagaactaataagagaaaaattccagctggagcaggagagaacaatctgctgaagacactggccagctcctgagcagcttgcaggtttcaggccagctgcagactggctccccctgccaccacagccacagggcaccacacagccccccagccagcagcagcagggctcagggctcaggaccaaccgcatag