

UQCR-ubiquinol-cytochrome c reductase, 6.4kDa subunit Gene

Reference: PTXBC000462

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UQCR; 0710008D09Rik; QCR10; cytochrome b-c1 complex subunit 10; complex III subunit 10; complex III subunit XI; ubiquinol-cytochrome c reductase complex 6.4 kDa protein; ubiquinol-cytochrome c reductase, 6.4kDa subunit; ubiquinol-cytochrome c reductase, c		

For research use only.

GENE DESCRIPTION	ubiquinol-cytochrome c reductase, 6.4kDa subunit
GENE ID	10975
ORIGIN SPECIES	Human
ACCESSIONS	BC000462
VECTOR	pDONR223
NCBI SYMBOL	UQCR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgacccggttcctgggccacgtaccgggagctggtcaagaactgggtcccgacggcctacacatggggcgctgtgggcgctgtggggcgtgggtgtgggccaccgattggcggtgatcctggactgggtaccttacatcaatggcaagtttaagaaggataattaa