

GADD45G-growth arrest and DNA-damage-inducible, gamma Gene

Reference: PTXBC000465

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CR6; DDIT2; GADD45gamma; GRP17; growth arrest and DNA damage-inducible protein GADD45 gamma; DDIT-2; DNA damage-inducible transcript 2 protein; GADD45-gamma; cytokine-responsive protein CR6; gadd-related protein, 17 kD; growth arrest and DNA damage induci		

For research use only.

GENE DESCRIPTION	growth arrest and DNA-damage-inducible, gamma
GENE ID	10912
ORIGIN SPECIES	Human
ACCESSIONS	BC000465
VECTOR	pDONR223
NCBI SYMBOL	GADD45G
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactctggaagaagtcgcggccaggacacagttccgaaagcacagccaggatgcagggtgccggaaagcgcgtgcatgagttgctgtcggcgacgctcagggtgctcactgcccggcgtctacgagtcagccaaagtctgaacgtggacccccgacaatgtgacctctgtgtgtcgtgctgcgggtgaggaggacgagggcgacatcgctgcagatccattttacgtgaltccaggctttctgctgcgagaacgacatcgacatagtcgctgtggcgatgtgcagcggctggcgctatctgtggcgccggcgagggagggcggtgcggcgacctgcac