

UBE2V1-ubiquitin-conjugating enzyme E2 variant 1 Gene

Reference: PTXBC000468

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CIR1; CROC-1; CROC1; UBE2V; UEV-1; UEV1; UEV1A; ubiquitin-conjugating enzyme E2 variant 1; DNA-binding protein; TRAF6-regulated IKK activator 1 beta Uev1A; ubiquitin conjugating enzyme E2 V1		

For research use only.

GENE DESCRIPTION	ubiquitin-conjugating enzyme E2 variant 1
GENE ID	7335
ORIGIN SPECIES	Human
ACCESSIONS	BC000468
VECTOR	pDONR223
NCBI SYMBOL	UBE2V1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagccaccacgggctcgggagtaaaatccctcgcaatttccgactgttggaagaactcgaagaaggccagaaggagtaggagatggcacagttagctgggtctagaagatgacgaagacatgacacttacaagatggacaggatgataattgggcctccaagaacaattta
tgaaccgaatatacagccttaaaatagaatgtggacctaataccagaagcacccttftgaagattgttaacaaaaataatgaatggagtaaatagttctaatggagtggtggaccaagagccatcatcagtgctagcaaaatggcagaattcatatagcatcaaaagtgtcctgca
agagcttcgggcctaatgatgtctaaagaaatagaaactccctcagccgccgaaggacagtggttacagcaattaa