

RPS4X-ribosomal protein S4, X-linked Gene

Reference: PTXBC000472

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCG2; DXS306; RPS4; SCAR; SCR10; 40S ribosomal protein S4, X isoform; cell cycle gene 2; ribosomal protein S4X isoform; single copy abundant mRNA protein; single-copy abundant mRNA; ribosomal protein S4, X-linked		

For research use only.

GENE DESCRIPTION	ribosomal protein S4, X-linked
GENE ID	6191
ORIGIN SPECIES	Human
ACCESSIONS	BC000472
VECTOR	pDONR223
NCBI SYMBOL	RPS4X
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctcgtgtcccaagaagcatctgaagcgggtggcagctccaaagcattggatgctgataaattgaccgggtgtgttctctctgctccatccaccgggtcccaagttgagagagtgctctcccctcatcttctgaggaacagacttaagtatgccctgacaggagatgaagtaa
agaagattgcatgcagcgggtcattaaaaatcgalggcaaggtccgaactgatataacctaccctgctggatcatggatgctatcagcattgacaagacgggagagaattccgctgctgatctgacaccaagggtcgttctgtacatcgtattacacfcaggaggccaagtacaagttg
tgcaaaftgagaaagatctttgtgggcacaaaaggaatccctcatctggtgactcatgatgcccgaccctccgclacccgatcccctcatcaagttgaatgataccattcagattgattggagactggcaagattactgattcattcaagttcgacactggtaacctgtgtatggtgactgga
ggtgctaactaggaagaattggtgatcaccacaagagagaggaccctggatctttgacgtggttcacgtgaaagatccaatggcaacagctttgccactgacttccaacattttgtattggcaaggccaacaacctgatttctctcccaggaaagggtatccgctcac
cattgctgaagagagagacaaaagactggcggccaacagagcagtggtgta