

RPL13A-ribosomal protein L13a Gene

Reference: PTXBC000514

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L13A; TSTA1; 60S ribosomal protein L13a; 23 kDa highly basic protein; tissue specific transplantation antigen 1; ribosomal protein L13a		

For research use only.

GENE DESCRIPTION	ribosomal protein L13a
GENE ID	23521
ORIGIN SPECIES	Human
ACCESSIONS	BC000514
VECTOR	pDONR223
NCBI SYMBOL	RPL13A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggaggtgcaggtcctggtgatggtcaggccatctcctggccgctggcgccatcgtgctaaacaggtactctgggcccgaaggtggtgctgtacgtgtgaaggcatcaattcttgcaatttctacagaacaagttgaagtacgtggttctcgcgaagc
ggatgaacaccaacccctcccgagggccctaccacttccggccccagccgcacattcttgccggaccgtgcgaggtatgctgccccacaaaacaaagcagggccagccgctctggaccgtctcaaggtgtttgacggcatccaccgcccctacgacaagaaaaagcggatggtg
ttcctgctgccctcaaggtcgtgctgaagcctacaagaaagttgctatctggggcgcctggctcacgaggttgctggaagtaccaggcagtgacagccacctggaggagaagaggaagagaagccaagatccactaccggaagaagaacagctcatgaggctacgga
aacaggccgagaagaacgtggagaagaaattgacaatacacagaggtcctcaagaccacggactcgtgctga