

RPS6-ribosomal protein S6 Gene

Reference: PTXBC000524

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	40S ribosomal protein S6; phosphoprotein NP33; ribosomal protein S6		
For research use only.			

GENE DESCRIPTION	ribosomal protein S6
GENE ID	6194
ORIGIN SPECIES	Human
ACCESSIONS	BC000524
VECTOR	pDONR223
NCBI SYMBOL	RPS6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagctgaacatctctccagccactggctgccagaaactcattgaagtggacgatgaacgcaactctgactttctatgagaagcgtatggccacagaagttgctgctgacgctctgggtgaagaatggaagggttatgtgtccgaatcagtggtgggaacgacaaacaaggttc
cccatgaagcagggtgtcttgacctggccgtgtccgcctgctactgagtaaggggcattcctgttacagaccaaggagaactggagaagaagagaaaaatcagttcgtgtgcatgttgatgcaaatctgagcgttctcaacttggtattgtaaaaaaggagagaaggatattcct
ggactgactgatactacagtgctcgcgcctgggccccaaaaagagctagcagaatccgcaaacctttcaatctctctaaagaagatgatgtccgcatgtgtgtaagaagcccttaataaagaaggtaagaacctaggaccaagcacccaagattcagcgtctgttactccacgt
gtcctgcagcacaacggcgcgtattgtctctgaagaagcagcgtaccaagaaaaataaagaagggtgcagaatatgctaaacttttggccaagagaatgaaggaggctaaaggagaagcgcagggaacaaattgcgaagagacgcagacttctctctgcgacttctacttctaa
tctgaatccagtcagaataa