



Reference: PTXBC000527

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EWS-FLI1; bK984G1.4; RNA-binding protein EWS; EWS RNA-binding protein variant 6; Ewing sarcoma breakpoint region 1; Ewings sarcoma EWS-Flil (type 1) oncogene; EWS RNA binding protein 1		

For research use only.

GENE DESCRIPTION	Ewing sarcoma breakpoint region 1
GENE ID	2130
ORIGIN SPECIES	Human
ACCESSIONS	BC000527
VECTOR	pDONR223
NCBI SYMBOL	EWSR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr