

ITM2B-integral membrane protein 2B Gene

Reference: PTXBC000554

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ABRI; BRI; BRI2; BRICD2B; E25B; E3-16; FBD; RDGCA; imBRI2; integral membrane protein 2B; ABri/ADan amyloid peptide; BRICHOS domain containing 2B; immature BRI2; transmembrane protein BRI		

For research use only.

GENE DESCRIPTION	integral membrane protein 2B
GENE ID	9445
ORIGIN SPECIES	Human
ACCESSIONS	BC000554
VECTOR	pDONR223
NCBI SYMBOL	ITM2B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgaaggtagcgttcaactccgctctggccagaaggaggccaagaaggagcgagcccaagagcggcgaggagcgctcatcatcccccgagccgctgcggtgggactgcaaggaccagatgatgtgtaccagttggccaaagaagagcctggtgttgcatgtcttt
ggactagcattatgcttcagggcgttattctaggagagcactgtgtacaaatatttgcactcaaccagatgacgtgtactactgtggaataaagtacatcaaaagatgatgtcatcttaafgagccctctgcagatccccagctgctctaccagacaattgaagaaatataaaatcttt
gaagaagaagaagttgaattatcagtgctgtccagagtttgcagatagtgatcctgccaacattgttcagactttaacaagaacttacagcctattagatcttaacctggataagtgctatgtgatccctctgaacacttccattgttatgccaccagaaacactaggagttactattaa
catcaaggctggaacctatttgcctcagtcctatctgattcatgagcacatggttattactatgcattgaaaacattgacacctgggtttctttatctgactgtgtcatgacaagaaactacaactgcaacgcagagaaactataaaggtattcagaacgtgaagccagcaattgttc
gcaattcggcattttgaaaacaaatttgcctgggaactttaattgttctga