

HIGD1A-HIG1 domain family, member 1A Gene

Reference: PTXBC000601

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HIG1; RCF1a; HIG1 domain family member 1A, mitochondrial; HIG1 domain family, member 1A; RCF1 homolog A; hypoxia inducible gene 1; hypoxia-inducible gene 1 protein; HIG1 hypoxia inducible domain family member 1A		

For research use only.

GENE DESCRIPTION	HIG1 domain family, member 1A
GENE ID	25994
ORIGIN SPECIES	Human
ACCESSIONS	BC000601
VECTOR	pDONR223
NCBI SYMBOL	HIGD1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcaacagacacaggtgttcccttccatcatgaggaagatcaggatcaaaactcattcgaaaagctaaaggacacattcgtaccgttggaatagcgggtttgcagcaattgtgcatatggattataaaactgaagacaggggaaatactaaaatgtccattcatctgatccac
atgcgtgtggcagcccaaggccttgttaggagcaatgactgttggtatggcctatccatgcatcgggaattctggcgaacctaagccttag