



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

RPL14-ribosomal protein L14 Gene

Reference: PTXBC000606

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

CAG-*ISL-7*; CTG-B33; L14; RL14; hRL14; 60S ribosomal protein L14; ribosomal protein L14

For research use only.

GENE DESCRIPTION	ribosomal protein L14
GENE ID	9045
ORIGIN SPECIES	Human
ACCESSIONS	BC000606
VECTOR	pDONR223
NCBI SYMBOL	RPL14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]