

ID1-inhibitor of DNA binding 1, dominant negative helix-loop-helix protein Gene

Reference: PTXBC000613

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	bHLHb24; DNA-binding protein inhibitor ID-1; class B basic helix-loop-helix protein 24; dJ857M17.1.2 (inhibitor of DNA binding 1, dominant negative helix-loop-helix protein); inhibitor of DNA binding 1, dominant negative helix-loop-helix protein; inhibito		

For research use only.

GENE DESCRIPTION	inhibitor of DNA binding 1, dominant negative helix-loop-helix protein
GENE ID	3397
ORIGIN SPECIES	Human
ACCESSIONS	BC000613
VECTOR	pDONR223
NCBI SYMBOL	ID1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaagtcgccagtggcagcaccgccaccgccgcgggccccagctgcgcgtgaaggccggcgaagacagcagcgggtgcggcgaggtggtgcgtgtctgtctgagcagagcgtggccatctcgcgtgcgcggggggcgccggggcgccgctgcctgcctgctgg
acgagcagcaggtaaacgtgctgctctacgacatgaacggctgttactacgcctcaaggagctggtgccccacctgccccagaaccgcaaggtgagcaaggtggagattctccagcagctcatcactacatcaggaccttcagttggagctgaactcggaatccgaagttggaacc
ccccggggcgccgagggctgccggctccggtcagcaccctcaacggcgagatcagcgcctgacggccgagggcgcatgcttctgcggacgatcgcatcttgtgtcgtga