

GRB2-growth factor receptor-bound protein 2 Gene

Reference: PTXBC000631

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	epidermal growth factor receptor-binding protein GRB2; SH2/SH3 adapter GRB2; EGFRBP-GRB2; ASH; Grb3-3; MST084; MSTP084; NCKAP2; growth factor receptor-bound protein 2; HT027; abundant SRC homology; growth factor receptor-bound protein 3; protein Ash; grow		

For research use only.

GENE DESCRIPTION	growth factor receptor-bound protein 2
GENE ID	2885
ORIGIN SPECIES	Human
ACCESSIONS	BC000631
VECTOR	pDONR223
NCBI SYMBOL	GRB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagccatcgccaaatgacttcaagctactgcagacgacgagctgagcttcaaaaggggggacatcctcaagggtttgaacgaagaatgtgatcagaactgttacaaggcagagcttaatggaaaagacggcttattccaagaactacatagaatgaaaccacatccgtggt
tttttgcataatccccagagccaaaggcagaagaatgcttagcaaacagcggcagatggggcctttctatccgagagagtgagagcgctcctggggacttccctctctgtaagtttggaaacgatgtgcagcacttcaagggtctccgagatggagccgggaagtacttctctgg
gtgtggaagtcaattcttgaatgagctgggtgattatcacagatctacatctgtctccagaaccagcagatattctctcgggacatagaacaggtgccacagcagccgacatacgtccagccctttgactttgattcccaaggaggtgagagctgggttccgcccgggagatttta
tccatgtcatggataactcagacccaactgggtgaaaggagcttgcacggggcagacggcatgtttcccgcaattatgtcaccctcgtaaccggaactctaa