

DHRS4L2-dehydrogenase/reductase (SDR family) member 4 like 2 Gene

Reference: PTXBC000663

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SDR25C3; dehydrogenase/reductase SDR family member 4-like 2; NADP(H)-dependent retinol dehydrogenase/reductase short isoform-like protein; NADPH-dependent retinol dehydrogenase/reductase-like protein 2; dehydrogenase/reductase (SDR family) member 4 like 2		

For research use only.

GENE DESCRIPTION	dehydrogenase/reductase (SDR family) member 4 like 2
GENE ID	317749
ORIGIN SPECIES	Human
ACCESSIONS	BC000663
VECTOR	pDONR223
NCBI SYMBOL	DHRS4L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaggctgctagcctctgtgctgggacggaagtcggtgctgggtggccagctccaggatgacccgccgggacccgctcacaataaggtggccctggtaacggcctccaccgacgggatcggtctcgccatgcccggcggtttggccaggacagggccacgtggtcgt
cagcagccggaagcagcagaatgtggaccaggcggtggccacgctgcagggggaggggctgagcgtgacgggacgtgtgccaatgtgggaaggcgaggacgggagcggtggtggccatggctgtgaagctcatggaggtatcgatctcctagctcctaatgctgctgic
aaccttctttggaagcctaattggtatcaccaggaggtgtgggacaagactctggacattaatgtgaaggccccagccctgatgacaaaggcagtggtgcccagaatggagaacgaggagggcggtcagtggtgatcgtgtcttccatagcagccttcagtcctcctggttca
gtccttacaatgtcagtaaaacagccttgcctggcctcaacaataccctggccatagagctggccccaaggacattagggtgaactgcctgcacctggacttatcaagactagcttcagcaggatgctctggatggacaaggaaaaaggaaagcatga