



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

CD320-CD320 molecule Gene

Reference: PTXBC000668

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD320 molecule; CD320 antigen; 8D6; 8D6A; TCBLR; 8D6 antigen; FDC-SM-8D6; FDC-signaling molecule 8D6; transcobalamin receptor		

For research use only.

GENE DESCRIPTION	CD320 molecule
GENE ID	51293
ORIGIN SPECIES	Human
ACCESSIONS	BC000668
VECTOR	pDONR223
NCBI SYMBOL	CD320
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacggcggtggtgatggcgacagtttgagcgtggcgacgaaggcgctctggcgctgctgctgctgctgcgcctcggactagtcctggagggccgcgagccgcttccacaccgcacctctccagagccgacggcccgacgtcgtgccaccaccacaagt
ccagtgccgcaccagtgtgcttatgctgctccctcacttgcgctgcgacagggactgtggactgcagcgtatggcagcgtatgaggagagtgtcaggattgagcactgtaccagaagaaggcgaatgccaccgcccccctggctcccttgccttgcaccggctcagtgtactgtcttgg
gggaactgacaaagaactctgcgaactcagcgcctggcctgctctagcagcgagctcctgtgacgtctgagcgtatgactgcatcctcactcactgcgcctgcgacggccacccactgtctgccagctccagcagcagctcgtgtgtgaacaatgagatctcccggaagggg
atgccacaacattggggccccctgtgacctgtgagatgtgctactctctcagaagtcgacacaaactctggggccccctgtgacctggagagctgtccctctctgggaatgccacactctctctgcggagaccctgtggaaggcccaactcctatggggttatgacgtctgcgtg
gtcgtctcaagcctgtgacacgcacccctctcttcttctctcgcagcgcagcgcctcggcctcgggtgtactctgtgctcattggcccatgaggaagtgctcctctgtcagacaagaagaactctcctctgta