



Reference: PTXBC000674

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DQ2; DYT1; torsin-1A; dystonia 1 protein; dystonia 1, torsion (autosomal dominant; torsin A); torsin A; torsin ATPase 1; torsin ATPase-1A; torsin family 1 member A		

For research use only.

GENE DESCRIPTION	torsin family 1, member A (torsin A)
GENE ID	1861
ORIGIN SPECIES	Human
ACCESSIONS	BC000674
VECTOR	pDONR223
NCBI SYMBOL	TOR1A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]