



Reference: PTXBC000683

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM46A; C6orf37; XTP11; HBV X-transactivated gene 11 protein; HBV XAg-transactivated protein 11; retinal expressed gene C6orf37; family with sequence similarity 46 member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 46, member A
GENE ID	55603
ORIGIN SPECIES	Human
ACCESSIONS	BC000683
VECTOR	pDONR223
NCBI SYMBOL	FAM46A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]