

RPL24-ribosomal protein L24 Gene

Reference: PTXBC000690

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEL-S-310; L24; 60S ribosomal protein L24; 60S ribosomal protein L30; epididymis secretory protein Li 310; ribosomal protein L30; ribosomal protein L24		

For research use only.

GENE DESCRIPTION	ribosomal protein L24
GENE ID	6152
ORIGIN SPECIES	Human
ACCESSIONS	BC000690
VECTOR	pDONR223
NCBI SYMBOL	RPL24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaggtcagctgtgcagtttttagcgggtacaagatctacccggacacgggagggcgtacgcccaggaccgacgggaaggtttccagtttcttaatgcgaaatgcgagtcgggtttctccaagggaatcctcggcagataaactggactgtctctacagaagggaagcaca
agggacagtcggaagaattcaaaagaaagaaacccgcccagcagtcacaaatccagagggccactactgtgcatctctctgatataatggccaagaggaatcagaacctgaagttagaaaggctcaacgagaacaagctatcagggtgctaaggaaagcaaaaggctaagca
agcatctaaaaagactgcaatggctgctgctaaggcacctacaaggcagcacctaagcaaaagattgtgaagcctgtgaaagtctcagctccccgagttgtgtgaaaacgctaa