

C21orf59-chromosome 21 open reading frame 59 Gene

Reference: PTXBC000709

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	c21orf59; kur; kurly; chromosome 21 open reading frame 59 L homeolog		
For research use only.			

GENE DESCRIPTION	chromosome 21 open reading frame 59
GENE ID	56683
ORIGIN SPECIES	Human
ACCESSIONS	BC000709
VECTOR	pDONR223
NCBI SYMBOL	C21orf59
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggttctgctgacgtgaagcggggcgacgagagccagttcctgctcaggcgctgggagtagcagctggaggagctcacggtgaggtggccgggtctataatggcggtcaagggtcagcgccctgctcagaaatggaagaattagccgaacatggcatatttcctcc
taatatgcaaggactgaccgatgatcagattgaagaattgaattgaaggatgaatgggtgaaaaatgcgtaccagcggaggtgagtggttaaaaaaggatgattggacgaaggatgggcaagctccaatgagaagatgaagcaagtgttaagaagactatagaagaagccaa
agcaataatctaaagaacaagtgaagccggtgctgtgttaccatggagatggtgaaagatgccttggaccagcttcgagcgccggtgatgattgtttaccccatgggtgtgccaccgtatgatccatccgcatggagtttgaaaaaaggaagacttgcgggaacacaggcagggc
tcaacgtcattaaagaggcagagcgcgagctgtggtggcgagccaaggagctgagaagaacgaagaagctttcagactacgtggggaagaatgaaaaaaccaaaattatcgcaagattcagcaaggggacagggagctccagcccgagagcctattattagcagtgaggagcag
aagcagctgatgctgtactatcacagaagacaagaggagctcaagagattggaagaaatgatgatgcctatttaactcaccatggcggtataacactgctttgaaaagacatttcatggagtgaagacataaaggagaccaagatga