



Reference: PTXBC000724

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	16A; ATP6IP1; ATP6S1; Ac45; CF2; VATPS1; XAP-3; XAP3; V-type proton ATPase subunit S1; ATPase, H+ transporting, lysosomal (vacuolar proton pump), subunit 1; ATPase, H+ transporting, lysosomal accessory protein 1; ATPase, H+ transporting, lysosomal interac		

For research use only.

GENE DESCRIPTION	ATPase, H ⁺ transporting, lysosomal accessory protein 1
GENE ID	537
ORIGIN SPECIES	Human
ACCESSIONS	BC000724
VECTOR	pDONR223
NCBI SYMBOL	ATP6AP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]