

CD70-CD70 molecule Gene

Reference: PTXBC000725

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD70 molecule; surface antigen CD70; CD70 antigen; CD27-L; CD27L; CD27LG; TNFSF7; TNLG8A; CD27 ligand; Ki-24 antigen; tumor necrosis factor ligand 8A; tumor necrosis factor ligand superfamily member 7		

For research use only.

GENE DESCRIPTION	CD70 molecule
GENE ID	970
ORIGIN SPECIES	Human
ACCESSIONS	BC000725
VECTOR	pDONR223
NCBI SYMBOL	CD70
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccggaggagggttcgggctgctcgtgcgcgcaggccctatgggtgcgtcctcgggctgcttggccattggtcgcgggcttggtgatctgcctcgtggtgtgcatccagcgccttcgcacaggctcagcagcagctgccgctcagtcacttgggtgggacgtagctgagctg
cagctgaatcacacaggacctcagcaggacccaggtatctgagcggggggccagcactggccgctccttcctgcatggaccagagctggacaaggggcagctacgtatccatcgtgatggcatctacatggtacacatccaggtgacgctggccatctgttctccacgacgg
cctccaggccaccaccaccaccctggcctgggaatctgctctccgctccgtagcatcagcctgctgcgtctcagcttcaccaagggtgtaccattgcctccagcgcctgacgcccctggcccaggggacacactctgcaccaacctcactgggacactttgcttcccga
cactgatgagaccttcttgagtgacgtgggtgcgcccctga