

CD82-CD82 molecule Gene

Reference: PTXBC000726

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD82 molecule; kangai 1 (suppression of tumorigenicity 6, prostate; CD82 antigen (R2 leukocyte antigen, antigen detected by monoclonal and antibody IA4)); CD82 antigen; 4F9; C33; GR15; IA4; KAI1; SAR2; ST6; TSPAN27; C33 antigen; inducible membrane protein		

For research use only.

GENE DESCRIPTION	CD82 molecule
GENE ID	3732
ORIGIN SPECIES	Human
ACCESSIONS	BC000726
VECTOR	pDONR223
NCBI SYMBOL	CD82
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtcagcctgtatcaaaagtcacaaatacttctctctctcaacttgatcttcttctctgggcgcagtgatctgggcttgggggtgtgatcctggccgacagagcagtttcatctctgtctgcaaacctctccagctcgttaggatggggcctatgtcttcacggcgtggg
ggcagtcactatgtcatgggcttctgggctgcatcggcgcgctcaacgaggtccgctgctgctggggtgtactttgcttctcgtctgactcctcattgccaggtgacggccggggccctcttctacttcaacatgggcaagctgaagcaggagatggggcgcacgtgactgagct
cattcgagactacaacagcagtcgcgagggacagcctgaggatgcttgggactacgtgcaggctcagggtgaagtgtcgcggctgggtcagcttcaactggacagacaacgtgagctcatgaatcgccctgaggtcacctacccctgttctcgcgaagtcaagggggaagaggac
aacagccttctgtgaggaagggttctcgcgagggcccggaacaggaccagagtggaaccaccctgaggactggcctgtgtaccaggagggtgcatggagaaggtgcaggcgtggctgcaggagaacctgggcatcatctcggcgtggcggtgtgtggccatcgtga
gctcctggggatggtcctgtccatctgtgtgcccgcacgtccattcgaagactacagcaaggtcccaagtactga