

EMD-emerin Gene

Reference: PTXBC000738

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EDMD; LEMD5; STA; LEM domain containing 5		
For research use only.			

GENE DESCRIPTION	emerin
GENE ID	2010
ORIGIN SPECIES	Human
ACCESSIONS	BC000738
VECTOR	pDONR223
NCBI SYMBOL	EMD
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacaactacgcagatcttccgataccgagctgaccacgtgctgcgcgggtacaacatccgcacggcgctgtagtaggatcaactcgtaggctttacgagaagaagatctcagtagcagaccagaggcgcggtctcgcggcggtcctcttatagct
tctctgactgaattcgactagaggggatgagatgatgatcttccaaagaaaggagcgtttactctacagagcaagggtacaatgacgactactatgaagagagctactcaccaccagagcttatggggagcccgagctcgcggcggtccagggtgctccgacgtagtg
acttcaatccagatgctgacgctttccatcaccaggtgcatgatgacgatctttgtctcttctgaagaggagtgcaaggatagggaacgcccagtgacggcgggacagtgccctaccagagcatcacgactaccgctgtttcagcctccaggagctccctggacctgtcctattatcc
tacttctctccacctctttatgtctctcctcatcttctctcatggctacccgcggtgccatccggcctgaaacccgtgctcctggggctgggctgggcccagatgccaggtcccgtctggggccagctgctgttttctggtctttgtgatgctctcttcttcaatccacttcatgc
aggctgaagaaggcaaccctcttag