

POLR2H-polymerase (RNA) II (DNA directed) polypeptide H Gene

Reference: PTXBC000739

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RPABC3; RPB17; RPB8; DNA-directed RNA polymerases I, II, and III subunit RPABC3; DNA-directed RNA polymerase II subunit H; DNA-directed RNA polymerases I, II, and III 17.1 kDa polypeptide; RPB8 homolog; polymerase (RNA) II (DNA directed) polypeptide H; po		

For research use only.

GENE DESCRIPTION	polymerase (RNA) II (DNA directed) polypeptide H
GENE ID	5437
ORIGIN SPECIES	Human
ACCESSIONS	BC000739
VECTOR	pDONR223
NCBI SYMBOL	POLR2H
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggcatcctgtttgaggatatttcgatgtgaaggatattgacccggagggaagaagtttgaccgagtgctcgcactgacattgtgagagtgaaatcttcaagatggatccttaattagatgtaaacattcaaatctaccctgtagacttgggtgacaagtgttcggttggtcatagctagtacc
ttgtatgaagatggtacccctggatgatgtggaatacaacccactgatgtaggcctccaggcctgaccagtttgagtafataatgtatggaagatgtacaggatgaggagatgaaacttctactgaagcagcaacacgcctctctgctacgtgtcctatggggcctgctcatgaggc
tgcaggggatgccacaacctgcatggattcgaggtggactccagagttatctctgatgaagaagctagcctctga