

RPL36AL-ribosomal protein L36a-like Gene

Reference: PTXBC000741

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RPL36A; 60S ribosomal protein L36a-like; ribosomal protein HL44; ribosomal protein L36a; ribosomal protein L36a like		

For research use only.

GENE DESCRIPTION	ribosomal protein L36a-like
GENE ID	6166
ORIGIN SPECIES	Human
ACCESSIONS	BC000741
VECTOR	pDONR223
NCBI SYMBOL	RPL36AL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtcaacgtacctaaccgaagaaccttctgtaagaagtgtggcaagcatcagctcacaaagtacacagtataagaaggcaaggattcttctgtatccccagggaaggaggcgctatgatcggaagcagagtggtctatgggtggcagacaaagccaatttccggaagaaggc
taagaccacaaagaagattgtgtaaggctggaatgtgttgagcctaactgcagatccaagaggatgctggccattaaagagatgcaagcatttgaactgggaggagataagaagagaaaggccaagtgtatccagttctaa