

EIF5A-eukaryotic translation initiation factor 5A Gene

Reference: PTXBC000751

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EIF-5A; EIF5A1; eIF5AI; eukaryotic translation initiation factor 5A-1; eIF-4D; eIF-5A-1; eIF-5A1; eukaryotic initiation factor 5A; rev-binding factor; eukaryotic translation initiation factor 5A		

For research use only.

GENE DESCRIPTION	eukaryotic translation initiation factor 5A
GENE ID	1984
ORIGIN SPECIES	Human
ACCESSIONS	BC000751
VECTOR	pDONR223
NCBI SYMBOL	EIF5A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagatgacttgactcgagacaggagatgcaggggcctcagccaccttcccaatgcagtgctcagcattacgtaagaatggctttgtggtgtctcaaggccgccatgtaagatcgtcgagatgtctactcgaagactggcaagcagggccacgccaaggtccatctgttggtat
tgacatcttactgggaagaatatgaagatatctgccgtcaactcataatgatggatgtccccaacatcaaaaggaaatgacttccagctgattggcatccaggatgggtacctatcactgtccaggacagcggggaggtacgagaggacctctctccctgaggagaccttggaagg
agattgagcagaagtacgactgtggagaagatcctgatcacgggtgctgtgcatgacagaggagggcagctgttgcaatcaaggccatggcaaaataa