

RAB13-RAB13, member RAS oncogene family Gene

Reference: PTXBC000799

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RAB13, member RAS oncogene family; RAS-associated protein RAB13; GIG4; ras-related protein Rab-13; cell growth-inhibiting gene 4 protein; growth-inhibiting gene 4 protein		

For research use only.

GENE DESCRIPTION	RAB13, member RAS oncogene family
GENE ID	5872
ORIGIN SPECIES	Human
ACCESSIONS	BC000799
VECTOR	pDONR223
NCBI SYMBOL	RAB13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaaagcctacgaccacctctcaagttgctgctgacgggggactcgggggtgggcaagacttgctgacattcgtttgcagaggacaactcaacaacacttacatctccaccatcggaattgattcaagatccgactgtggatatagaggggaagaagatcaaaactacaagctctgggacacggctggccaagagcgggttcaagacaataactactgcctactacgtggagccatgggcatatcctagtatacgacatcacggatgagaaatcttcgagaataticagaactggatgaaaagcatcaaggagaatgcctcggctgggtggagcgctcttctggggaacaatgtgacatggaggccaagaggaaaggtgcagaaggagcagggcgataaagttggctcgagagcatggaatccgattttcgaactagtgctaaatccagatgaatgtggatgaggcttttagttccctggccgggacatctgctcaagtcaggagccggagatcaggaaacggcaacaagcctcccagtagtgcctgaaaactgtgacaagaagaacacacaagtgctcctgggctga