



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SLC30A5-solute carrier family 30 (zinc transporter), member 5 Gene

Reference: PTXBC000808

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZNT5; ZNTL1; ZTL1; ZnT-5; zinc transporter 5; solute carrier family 30 (zinc transporter), member 5; zinc transporter ZTL1; znT-like transporter 1; solute carrier family 30 member 5		

For research use only.

GENE DESCRIPTION	solute carrier family 30 (zinc transporter), member 5
GENE ID	64924
ORIGIN SPECIES	Human
ACCESSIONS	BC000808
VECTOR	pDONR223
NCBI SYMBOL	SLC30A5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggagaataatcgcgggggacgtctgcggcgcccgccggcgggcgccgttggcgcggtggagctaccacgcgtcgtgattacaacaaatattgtgtactatgttctactaaattttgaaggctgtgggactttcgaatcatatgatctctaaagactgttcacattgttcagtt
ctatttatataaaacttggactgcatttatgtgttttttcaaaaccattttctctgcgaacaaattaccacaaccagataatttgatcactaaaaattcctgtgaaaaaaatttgaagacaaaaattgaatgatcctaggaactatgtgcgaactga