

FBXW5-F-box and WD repeat domain containing 5 Gene

Reference: PTXBC000850

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Fbw5; F-box/WD repeat-containing protein 5; F-box and WD-40 domain-containing protein 5; WD repeat-containing F-box protein FBW5; F-box and WD repeat domain containing 5		

For research use only.

GENE DESCRIPTION	F-box and WD repeat domain containing 5
GENE ID	54461
ORIGIN SPECIES	Human
ACCESSIONS	BC000850
VECTOR	pDONR223
NCBI SYMBOL	FBXW5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggcctgtcggcgacaacaggtacgtgtacgtgaacagccgcctggcccaacgggtgcgtgtgtggcgaccccatgcagccgcaccaatcgaggagattgacgtgctgtgttcgacctcaagaccatcggggaggtgaggcgggctctgctgcgcaccgcgcta
cagcccaacgacgagtgcttctcatcttctggagtcagcagggaactctgtggccagcggggcgaggaccgcacggctacatctgggaccgcaactacaacatctgtctggccaggtcgggcagcaggatgtgtcaactcagtggtcttcagtcgccaggagcaggagctg
ctgtcacggccagcgacgacgccaccatcaaagcctggcgtctccacgacacatgcgcgtctccaggcacctcgccacggcctgcacctcttctctgtgttgcagcagagggcgtga