

SSNA1-Sjogren syndrome nuclear autoantigen 1 Gene

Reference: PTXBC000864

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	N14; NA-14; Sjogren syndrome nuclear autoantigen 1; Sjogren syndrome nuclear autoantigen 1; Sjogren's syndrome nuclear autoantigen 1; nuclear autoantigen of 14 kDa; SS nuclear autoantigen 1		

For research use only.

GENE DESCRIPTION	Sjogren syndrome nuclear autoantigen 1
GENE ID	8636
ORIGIN SPECIES	Human
ACCESSIONS	BC000864
VECTOR	pDONR223
NCBI SYMBOL	SSNA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccacagcagggcgccgctgcagaactacaacaacagctgtgtaagtcataagaggagctgtgccagaagcgggaggagctgtgccggcagatccaggaggaggagcagaagcagcggctgcagaatgaggtgaggcagctgacagagaagctggcccgctca
acgagaacctggcagcgaagattgcctctcgcaacgagttcgaccggaccatcgccggagacggagcgcctacctcaagatcctggagagctcccagactttgctcagcgttctcaagagggaagctgggaacctgaccaaggctacagccccagaccagaaaaagtagcggcggc
agggacagctga