

MGC21881-hypothetical locus MGC21881 Gene

Reference: PTXBC000869

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized locus MGC21881		
For research use only.			

GENE DESCRIPTION	hypothetical locus MGC21881
GENE ID	389741
ORIGIN SPECIES	Human
ACCESSIONS	BC000869
VECTOR	pDONR223
NCBI SYMBOL	MGC21881
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagaaagctgagactcagggccagcaacccgggtcccagcggagcggccggcacacggcgcacttcagcaccagtggcgggtggccaccactgtgcgggagatggctgcgacgcgtgcgcagatctcgatccaaactccctcctgccagaatctggaccggaatccaccat
tgctcgtttccgctgcgcgtggagagaatctctgaggtccccaggagagcctgcctgcacggaagagacgcctcctcgattaaatgctgtttccatga